


Tasterbausatz M3 XXT, groß

626103-9025-000

Pos.	Art. Nr.	Stk.	Bezeichnung	Gewicht
1	626103-6110-004	1	Klemmschraube M3 XXT, L8, DG3,5	0,5 g
2	626103-6110-005	2	Klemmschraube mit Innengewinde M3 XXT, L8, DG3,5	0,7 g
3	626109-0026-000	1	XXT Taster Montagehilfe, Aluminium, Außenseite gerändelt	1,9 g
4	000000-0015-323	3	Stiftschlüssel für Tastersysteme M3 XXT, 1,1, CZ-Norm 1170	0,6 g
5	626103-6110-001	1	Würfel M3 XXT, 90 Grad, einfach, L6,5, B5	0,6 g
6	626103-6110-002	1	Würfel M3 XXT, 180 Grad, zweifach, L13, B5	1,1 g
7	626103-6110-003	1	Würfel M3 XXT, Stern, L13, B5	1,8 g
8	626103-6180-001	1	Gelenkstück M3 XXT, Titan, B5, L10	1,7 g
9	626103-5100-000	1	Tasterverlängerung M3 XXT, ThermoFit, DGE5, LE100	2,1 g
10	626103-5050-000	1	Tasterverlängerung M3 XXT, ThermoFit, DGE5, L50	1,5 g
11	626103-0500-075	1	Taster M3 XXT, ThermoFit Schaft, Siliziumnitridkugel, DG5, DK5, DS3,5, L75, ML66	1,8 g
12	626103-0300-050	1	Taster M3 XXT, ThermoFit Schaft, Rubinkugel, DG5, DK3, DS2, L50, ML41	0,8 g
13	626103-0300-033	2	Taster M3 XXT, ThermoFit Schaft, Rubinkugel, DG5, DK3, DS2, L33, ML24	0,7 g
14	626103-0500-030	2	Taster M3 XXT, ThermoFit Schaft, Rubinkugel, DG5, DK5, DS3,5, L30, ML21	1,1 g
15	626113-0200-020	2	Taster M3 XXT, Hartmetall Schaft, Rubinkugel, DG5, DK2, DS1,0, L20, ML11	0,7 g
16	626113-0150-020	2	Taster M3 XXT, Hartmetall Schaft, Rubinkugel, DG5, DK1,5, L20, ML11, DS1	0,7 g
17	626113-0100-020	2	Taster M3 XXT, abgesetzt, Hartmetallschaft, Rubinkugel, DG5, DK1,0, L20, ML11, MLE5, DS1, DSE0,8	0,7 g



Stylus kit M3 XXT, large

626103-9025-000

Pos.	Item No.	Qty.	Description	Weight
1	626103-6110-004	1	Clamp screw M3 XXT, L8, DG3,5	0,5 g
2	626103-6110-005	2	Clamp screw with female thread, M3 XXT, L8, DG3,5	0,7 g
3	626109-0026-000	1	XXT styli assembling aid, aluminum, knurled	1,9 g
4	000000-0015-323	3	Pin wrench, 1,1, CZ-Norm 1170	0,6 g
5	626103-6110-001	1	Cube M3 XXT, 90° Grad, single, L6,5, B5	0,6 g
6	626103-6110-002	1	Cube M3 XXT, 180°, double L13, B5	1,1 g
7	626103-6110-003	1	Cube M3 XXT, star L13, B5	1,8 g
8	626103-6180-001	1	Knuckle joint M3 XXT, titanium, B5, L10	1,7 g
9	626103-5100-000	1	Stylus extension M3 XXT, ThermoFit, DGE5, LE100	2,1 g
10	626103-5050-000	1	Stylus extension M3 XXT, ThermoFit, DGE5, L50	1,5 g
11	626103-0500-075	1	Stylus M3 XXT, ThermoFit shaft, silicon nitride sphere, DG5, DK5, DS3,5, L75, ML66	1,8 g
12	626103-0300-050	1	Stylus M3 XXT, ThermoFit shaft, ruby sphere, DG5, DK3, DS2, L50, ML41	0,8 g
13	626103-0300-033	2	Stylus M3 XXT, ThermoFit shaft, ruby sphere, DG5, DK3, DS2, L33, ML24	0,7 g
14	626103-0500-030	2	Stylus M3 XXT, ThermoFit shaft, ruby sphere, DG5, DK5, DS3,5, L30, ML21	1,1 g
15	626113-0200-020	2	Stylus M3 XXT, tungsten carbide shaft, ruby sphere, DG5, DK2, DS1,0, L20, ML11	0,7 g
16	626113-0150-020	2	Stylus M3 XXT, tungsten carbide shaft, ruby sphere, DG5, DK1,5, L20, ML11, DS1	0,7 g
17	626113-0100-020	2	Stylus M3 XXT, stepped, tungsten carbide shaft, ruby sphere, DG5, DK1,0, L20, ML11, MLE5, DS1, DSE0,8	0,7 g